

Lentiviral Production and Transduction Procedure

I. Plate HEK cells

- a. Coat plate with Poly-L-Lysine
 - i. Remove excess and allow plate to dry
- b. Seed plate with HEK cells (in D10)
 - i. 18×10^6 cells / 15cm dia. Dish
 - ii. Incubate overnight, 37C
 - iii. Plate should be 70%<x<90% confluent before transfection

II. Transfection

- a. Prepare DNA in serum-free media
 - i. 20ug/15cm plate: **packaging signal** (PCDNLBH) and 20ug/15cm **vector DNA**
 - ii. 5ug/15cm plate: **packaging envelope** pMN-VSV-G (or 1.5-2ug/15cm plate EO1-Foamy) for lentivirus
- b. Prepare Lipofectamine in serum-free media, incubate 5min @ RT
 - i. Add Lipofectamine dropwise and swirl gently
 1. 80ul/15cm plate
- c. Combine DNA and Lipofectamine solutions and incubate 30min @ RT
- d. Remove media from HEK cells (just before end of incubation)
- e. After incubation period, add mixed solution to HEK cells + 10ml serum-free media for 1-4hrs @ 37C
- f. Add 15ml growth media and incubate overnight @37C
- g. Change to growth media with 20ul .5M sodium butyrate/ml media, incubate 6-8hrs
- h. Change to growth media, incubate 18-24hrs.

D10 (growth media)	DMEM, 10% FBS, 1% Pen/Strep, 1% Sodium Pyruvate, 1% Non Essential Amino Acids
------------------------------	---

III. Virus Harvest

- a. Lenti Virus:
 - i. Remove and filter media from HEK cell plates
 1. Spin down media prior to filtering to pellet dead cell debris
 - a. 400 RCF
 - b. 5 min
 2. Use .45 vacuum filter
 - ii. Using *Centricons* spin down virus media and collect in filter
 1. Max 60 ml in *Centricons*
 2. 3000g for 30 minutes
 3. Repeat with any other media
 - iii. Remove virus from filter
 1. remove *Centricons* filter apparatus
 2. turn upside down

3. attach to conical collection dish
4. re-spin 1000g for 15min(or 10min)
- iv. Collect virus
 1. Remove virus from conical collection tube
 2. Determine M.O.I.
- v. Freeze @ -80°C
- b. Retro Virus:
 - i. Remove virus media from HEK cell plates.
 - ii. Centrifuge to remove cell debris
 - iii. Filter supernatant with .45 vacuum filter
 - iv. Determine M.O.I.
 - v. Flash freeze with Liquid N₂
 - vi. Store @ -80°C

IV. **Retronectin Transduction**

- a. Dissolve Retronectin lyophilized powder in sterile water to a concentration of 1mg/ml by gentle swirling (do not vortex) and filtrate through 0.22µm filter (Millipore Millex GV).
- b. Store solution in aliquots at -20°C (can be stored for one year).
- c. Coating:
 - i. Use only non-tissue culture-treated six-well plates (Becton Dickinson 351146, polystyrene) otherwise it won't work.
- d. Dilute the 1mg/ml stock solution of Retronectin to a concentration of 25-50µg/ml in sterile PBS.
- e. Dispense 2ml of the diluted Retronectin into each well and incubate the covered plate for 2 hours at room temperature (or at 4°C overnight).
- f. Remove Retronectin solution and block each plate with 2 ml PBS (containing 2% BSA) for 30 minutes at RT.
- g. Remove PBS (+2% BSA) and wash the wells once with 2 ml PBS. (The coated plate can be stored at 4 °C for one week)
- h. Remove the PBS and add virus supernatant (add growth media if needed to just cover surface of well)
- i. Incubate virus on retronectin for 15 min
- j. Add 1×10^6 Cells /cm² in each well with growth medium and centrifuge for 30 min with 600xg at 32°C.
- k. Incubate the cell for 2-days in the incubator (32-37°C).